uniGAV[®]

Ⓞ Instructions for use | Ⓞ 使用适应症

Ⓞ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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INDICATION	4
TECHNICAL DESCRIPTION	4
PHYSICS BACKGROUND	5
FUNCTION OF THE VALVE	6
POSSIBLE SHUNT COMPONENTS	6
TUBE SYSTEMS	7
SURGICAL PROCEDURE	7
PREOPERATIVE VALVE TEST	8
PRESSURE-FLOW CHARACTERISTICS	8
RE-IMPLANTATION	8
SAFETY MEASURES	8
COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	8
POSTOPERATIVE VALVE TEST	9
FUNCTIONAL SAFETY	9
STERILISATION	9
RESTERILISATION	9
NOTE ON THE INSTRUCTIONS FOR USE	9
MEDICAL PRODUCTS CONSULTANT	9
REQUIREMENTS OF THE MDD 93/42/EEC	10
NOTE ON THE INSTRUCTIONS FOR USE	10
GENERAL INFORMATION	10
VARIATIONS	11

INDICATION

The uniGAV is used for draining cerebrospinal fluid (CSF) from the ventricles into the peritoneum in hydrocephalus patients.

TECHNICAL DESCRIPTION

The uniGAV is a posture dependent hydrocephalus valve. It consists of a ball-cone unit and a gravitational unit. This ensures that optimal CSF drainage is achieved for each individual patient in any body position.

Fig.1 shows a schematic cross section of the uniGAV.

The uniGAV is made of a solid titanium body with a well-ried ball-cone unit integrated in its proximal part. A spiral spring (3) defines the opening pressure of the ball-cone unit.

The gravitational unit in the distal part contains a tantalum ball (5), which defines the opening pressure of this valve, and a sapphire ball (6), which ensures the precise closure of the valve.


Fig. 1: Schematic cross section of the uniGAV

PHYSICS BACKGROUND

The intraventricular pressure is positive in a healthy human in a horizontal position. To adjust this pressure through shunt drainage, one has to choose the appropriate pressure range, taking into account the abdominal cavity pressure. The resulting IVP is the sum of the shunt opening pressure and the abdominal cavity pressure (fig. 2).

The ventricular pressure in a healthy human in a vertical position becomes slightly negative. To maintain this pressure by means of shunt drainage, the shunt opening pressure has to be significantly higher so that the shunt can compensate the hydrostatic pressure minus the sum of the abdominal cavity pressure. Conventional shunts open immediately as soon as the patient stands up, which can lead to critical overdrainage.


Fig. 2: Calculating the intraventricular pressure

FUNCTION OF THE VALVE

The opening pressure of the *uniGAV* is composed of the opening pressure of the ball-cone unit and the opening pressure of the gravitational unit.

Horizontal position

The operational principle of the *uniGAV* is illustrated in figures 3a and b.

When the patient is in lying position, the gravitational unit is always open and therefore does not present any resistance to the fluid flow.

In fig. 3a, the ball-cone valve is closed. The drainage is blocked. In fig. 3b, the ball-cone unit is shown in the open position. The patient's IVP is increased and the spring force, which otherwise keeps the ball-cone unit closed, is overcome. The closing ball moves out of the cone and a gap opens to allow drainage. The opening pressure of the ball-cone valve is 7 cmH₂O.


Fig. 3: ball-cone unit in horizontal position
a) closed b) open

Vertical position

The operation of the gravitational valve is position-dependent. As soon as the patient moves into an upright position, the gravitational valve closes, the opening pressure of the *uniGAV* is significantly increased and CSF drainage is blocked (fig. 4a). Only when the sum of the IVP and the hydrostatic pressure exceeds the opening pressure of the *uniGAV*, drainage will be possible again (fig. 4b). The opening pressure of the *uniGAV* in the vertical position is 27 cmH₂O and the sum of the opening pressures of both the ball-cone valve and the gravitational valve (weight of the tantalum ball).

The total opening pressure refers to a reference flow of 5 ml/h. When the flow reach 20 ml/h, the opening pressures are approximately 1-2 cmH₂O higher.


Fig. 4: Gravitational unit in vertical position
a) closed b) open

POSSIBLE SHUNT COMPONENTS

The *uniGAV* is available with different shunt accessories. These variants are comprised of a variety of components, which are described below. There are variations for the use by children and adult hydrocephalus patients.

The *borehole reservoir* is positioned in the cranial borehole. It allows measuring the intraventricular pressure, injecting drugs and extracting CSF. Its solid titanium base is highly puncture-resistant. All reservoirs are available with integrated catheters or connectors. The *SPRUNG RESERVOIR* is a special borehole reservoir. As additional new feature of this reservoir CSF can be flushed towards the valve because of a one-way valve in the bottom of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR*.

The *prechamber* is positioned on the cranium. It allows measuring the intraventricular pressure, injecting drugs, extracting CSF and performing a palpatory ventricle inspection. Its solid titanium base is highly puncture-resistant. A puncture of the prechamber or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of max. 0,9 mm. 30 times of

punctures are able without any restrictions. A special prechamber is the *CONTROL RESERVOIR*. As an additional new feature of this reservoir, CSF can be flushed towards the valve because of a one-way valve in the proximal inlet of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.

Warning note: Frequent pumping can lead to overdrainage and thus to unphysiological pressure conditions. The patient should be informed about the risk.

Due to its tight fit on the ventricular catheter, the *deflector* allows choosing the length of catheter penetrating into the skull prior to implantation. The ventricular catheter is deflected at a right angle in the borehole (see chapter "Variations").

TUBE SYSTEMS

The *uniGAV* has been designed to ensure the optimal ventricular pressure. It is available as a shunt system or as individual valve units with or without an integrated distal catheter (internal diameter 1.2 mm, external diameter 2.5 mm). Individual valve units should be used with catheters of approx. 1.2 mm internal diameter and approx. 2.5 mm external diameter. The connector on the valve allows using catheters of 1.0 mm to 1.5 mm internal diameter. The external diameter of the catheter should be about double the internal diameter. In any case, the catheters must be carefully fixed, with a ligature, to the valve connectors. Kinks in the catheter have to be avoided.

The provided catheters have virtually no effect on the Pressure-flow characteristics.

SURGICAL PROCEDURE

Positioning the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The necessary skin incision should be carried out, preferably, in the shape of a lobule pedicled towards the draining catheter. To avoid CSF leakage, care should be taken that the dura opening is kept

as small as possible after applying the borehole. The ventricular catheter is stiffened by the introducing stylet supplied with the product.

The *uniGAV* is available in different shunt variants:

When using a *uniGAV SHUNTSYSTEM with borehole reservoir* or *SPRUNG RESERVOIR*, the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if CSF is dripping out. The catheter is shortened and the borehole reservoir is connected, with the connection secured with a ligature. The skin incision should not be located directly above the reservoir.

The *uniGAV SHUNTSYSTEM with prechamber* or *CONTROL RESERVOIR* comes with a deflector. This *deflector* is used for adjusting the position of deflection before implantation of the ventricular catheter. The catheter is deflected; the prechamber is put into place. The position of the ventricular catheter should be inspected again by postoperative CT or MR imaging.

Positioning the valve

The *uniGAV* is a posture-dependent shunt. Therefore, care must be taken that the valve is implanted parallel to the body axis. A suitable implantation site is behind the ear. The site of the implantation has no influence to the function of the valve.

After the skin incision and tunneling under the skin, the catheter is pushed forward, from the borehole to the intended shunt implantation site. The catheter is shortened, if necessary, and secured at the *uniGAV* with a ligature. The shunt should not be located directly under the skin incision.

The valve is marked with an arrow pointing in the direction of flow (arrow pointing to distal or downward).

Caution: We recommend for the connection of a catheter the use of armored clamps. Frequent pumping of a reservoir can lead to overdrainage and thus to unphysiological pressure conditions. The patient should be informed about the risk.

Positioning the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. It can be applied e. g. para-umbilically in a horizontal direction

or transrectally at the height of the epigastrium. Various surgical techniques are available for positioning the peritoneal catheter.

We recommend pulling through the peritoneal catheter, by means of a subcutaneous tunneling tool and perhaps with an auxiliary incision, from the shunt to intended position of the catheter. The peritoneal catheter, which is usually securely attached to the *uniGAV*, has an open distal end, but no wall slits. Following the exposure of, and the entry into, the peritoneum using a trocar, the peritoneal catheter (shortened, if necessary) is pushed forward into the open space in the abdominal cavity.

PREOPERATIVE VALVE TEST


Fig. 5: Patency test

The *uniGAV* can be filled by aspiration through a sterile, single-use syringe attached to the distal end of the catheter. The proximal end of the valve is immersed in a sterile, physiological saline solution. The valve is patent if fluid can be extracted in this way (see Fig. 5).

Caution: Pressure admission through the single-use syringe should be avoided, both at the proximal and the distal end.

Contaminations in the solution used for the test can impair the product's performance.

PRESSURE-FLOW CHARACTERISTICS

The following diagrams show the pressure-flow characteristics of the adjustable valve of the *uniGAV* for the pressure settings 7/27 cmH₂O. The total opening pressure refers to a reference flow of 5 ml/h. When the flowrates reach 20 ml/h, the opening pressures are approximately 1-2 cmH₂O higher.


Fig. 6: Pressure-flow characteristics for the pressure settings of the *uniGAV*

RE-IMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

SAFETY MEASURES

The patients must be carefully monitored after the implantation. Reddening of the skin and inflammation in the area of the drainage tissue could indicate infections at the shunt system. Symptoms such as headache, dizzy spells, mental confusion or vomiting are common occurrences in cases of shunt dysfunction. Such symptoms, as well as shunt system leakage, necessitate the immediate replacement of the shunt component responsible or of the entire shunt system.

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried out without endangering or impairing the functionality of the shunt.

The *uniGAV* is MR Conditional (ASTM F2503-08). All components are visible via X-ray. The provided catheters are MRI Safe. Reservoirs, deflectors and connectors are MR Conditional.

POSTOPERATIVE VALVE TEST

The *uniGAV* has been designed as a safe and reliable unit even without the provision of a pumping device. However, there are ways of testing the valve if a shunt system with a *prechamber* or a *borehole reservoir* is used. Valve tests can be carried out by flushing or pressure measurements.

FUNCTIONAL SAFETY

The valves have been designed for long-term reliable and precise operation. Still, the possibility that the shunt system needs to be replaced for technical or medical reasons, cannot be excluded. The valve and the valve system are able to resist positive and negative pressure up to 200 cmH₂O during and after implantation.

STERILISATION

The products are sterilised with steam under closely monitored conditions. The double wrapping in sterile bags ensures sterility for a period of five years. The expiry date is printed on the wrapping of each individual product. Products taken from a damaged wrapping must not be used under any circumstances.

RESTERILISATION

The functional safety and reliability of resterilized products cannot be guaranteed, therefore re-sterilisation is not recommended.

NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

MEDICAL PRODUCTS CONSULTANT

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH&Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz

Christoph Miethke GmbH & Co. KG
Ulanenweg 2
D-14469 Potsdam
Phone: +49(0) 7000 6438453 or
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
e-mail: info@miethke.com

Please address any enquiries to:

AESCULAP AG
Am Aesculap Platz
D-78532 Tuttlingen
Phone: +49 (0) 7461 95 9
Fax: +49 (0) 7461 95 26 00
e-mail: info@aesculap.de

Service address in the US
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042

AESCULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada
AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1-800-282-9000
www.aesculapusa.com

REQUIREMENTS OF THE MDD 93/42/EEC

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

GENERAL INFORMATION

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	<i>uniGAV</i>
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Drawing of the <i>uniGAV</i> with its external dimensions:	
Maßstab: 1:1	

VARIATIONS

The uniGAV is available as a single valve or as a shunt system comprising various components.

uniGAV

uniGAV SHUNTSYSTEM

uniGAV SHUNTSYSTEM with SPRUNG RESERVOIR or borehole reservoir (adult and pediatric)

uniGAV SHUNTSYSTEM with CONTROL RESERVOIR or prechamber (adult and pediatric)

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany